


Protein Electrophoresis & Transblot System

FL-ZY03



- FL-ZY03 Transfer Electrophoresis Tank (Holds 2 Gels)

Features

- Rapid, high-quality mini protein gel blotting
- Blots two 8.3×7.3cm gels at once
- The tank is integrally injection-molded from high-transparency polycarbonate material, ensuring leakproofness and impact resistance.
- 4cm electrode spacing: strong electric field for efficient transfer
- 99.99% Pt electrodes: better conductivity, corrosion/aging resistant
- Built-in blue crystal ice box: fast heat dissipation
- Open-type sandwich clamps: easy operation
- Includes bubble roller for de-bubbling
- Transfer area: 110×90mm
- Dimensions: 16×12×18cm
- Buffer volume: 1100mL

Standard Configuration

Name	Unit	Quantity
Electrophoresis Tank Lid (with cables)	SET	1
Electrophoresis Tank Lower Casing	PC	1
Transfer Inner Core	SET	1
Sandwich Transfer Clamp	SET	2
Transfer Sponge Pad	PC	4
Blue Ice Box	PC	1
Bubble Roller	PC	1